

Thermo Scientific Solaris qPCR Gene Expression Master Mix Plus ROX Vial

Introduction

Solaris qPCR Master Mix has been developed in conjunction with Solaris qPCR Gene Expression Assays for quantification of cDNA. Optimal qPCR results using probe detection chemistries can be achieved by using this system of qPCR products. These optimized master mixes offer a unique feature – they're blue. Working with blue reagents allows you to see your reagent as you pipette, track your progress across multiple wells, and visually assess your set-up. This visual confirmation further enhances the repeatability of your data. Solaris master mixes are available for any real-time instrument platform.

With the exception of probe, primers and template the Solaris qPCR Gene Expression Master Mix contains all the components required for reliable, sensitive and specific qPCR detection. These components include:

- Thermo-Start™ DNA Polymerase – a chemically modified hot-start version of ThermoPrime DNA Polymerase, which prevents non-specific amplification during the reaction set-up. **This enzyme performs optimally with a 95°C for 15 min activation step for complete enzyme activation.**
- Proprietary reaction buffer – a highly optimized buffering system that delivers improved assay performance when used in conjunction with Solaris primer/probe assays. The Solaris Master Mix also contains an inert blue dye to significantly increase the contrast between the reagent and plastic, making verification of master mix dispensing quick, easy and foolproof.
- dNTP – dTTP completely replaces dUTP in the Solaris qPCR Gene Expression Master Mix for maximum amplification efficiency.
- ROX – a separate vial of the internal passive reference dye ROX is supplied with this master mix and should be added at the required concentration specific to each qPCR instrument to allow normalization of data on ROX-dependent qPCR platforms.

Kit Contents

Vial	Pack Size (cap color)			
	INT (100 x 25 µL rxns)	A (200 x 25 µL rxns)	B (400 x 25 µL rxns)	C (1000 x 25 µL rxns)
Solaris qPCR Master Mix (2X)	1 x 1.25 mL (clear)	2 x 1.25 mL (clear)	4 x 1.25 mL (clear)	2 x 6.25 mL (clear)
ROX Reference Dye (1 mM)	25 µL (brown)	25 µL (brown)	25 µL (brown)	25 µL (brown)

Information

ROX Dye

ROX is an internal passive reference dye used to normalize the fluorescent reporter signal generated in qPCR. A separate vial of ROX is included for optional addition to the Solaris qPCR Gene Expression Master Mix. The final concentration will vary depending on the real-time instrument. For example, for a concentration of 100 nM ROX in a final 1X qPCR reaction mix, dilute ROX (1 mM) 40 times i.e. 5 µL ROX Reference Dye + 195 µL PCR grade water and add 10 µL of the diluted ROX solution to each 1.25 mL vial of Solaris Master Mix or 40 µL to each 5 mL vial of Solaris qPCR Master Mix.

Cycler Compatibility

Solaris qPCR Gene Expression Master Mix Plus ROX Vial is compatible with the qPCR instruments from all major suppliers.

Template Preparation

For maximum sample purity, we recommend extracting total RNA using a spin column method. We also

recommend using the Thermo Scientific Verso cDNA Synthesis Kit in the upstream reverse transcription (RT) step, as we have verified that this kit performs well with Solaris reagents.

Storage Conditions

Store at -20°C until ready for use. Solaris qPCR Gene Expression Low ROX Master Mix is stable for a minimum of 12 months. The reagents can be stored at 4°C for up to 1 month. Avoid repeated freeze thawing. The ROX dye is light sensitive; exposure should be minimized. Shipped on ice internationally and within the US.

Solaris qPCR Gene Expression Assays

Solaris qPCR Gene Expression Master Mixes are designed to be used with predesigned Solaris qPCR Gene Expression Assays for reliable, sensitive and specific detection.

Directions for Use

qPCR Protocol

1. Thaw the reagents on ice, mix the solutions and spin down before use to recover the maximum amount. **Do not vortex the Solaris Master Mix.**
2. Reaction mix preparation:

		96-well Plates	384-well Plates
		Volume	Volume
Reaction Mix	Solaris qPCR Master Mix (2X)	12.5 µL	5 µL
	Solaris Primer/Probe Set (20X)	1.25 µL	0.5 µL
	Water (PCR grade) ¹	Variable	Variable
	cDNA Template (1-250 ng) ²	1 - 5 µL	1 - 2 µL
Total volume		25 µL	10 µL

3. Seal plate with optically clear seal and briefly centrifuge to avoid bubbles within the wells, as these will interfere with the fluorescence detection. Always include a no template control (NTC) to assess any potential reagent contamination.
4. Run on thermal cycler.

qPCR Thermal Cycling Program:

	Temp.	Time	Number of cycles
Enzyme activation	95°C	15 min ³	1 cycle
Denaturation	95°C	15 sec	40 cycles
Annealing/Extension ⁴	60°C	60 sec	

Notes

1. The volume of the total reaction should be completed up to 10 µL or 25 µL with PCR grade water.
2. The volume of template added can be adjusted as required. For standard templates only 1 µL should be added to reduce carryover of PCR inhibitors. This volume can be increased for low copy number templates.
3. This enzyme performs optimally with a 95°C for 15 min activation step for complete enzyme activation.
4. Fluorescent data to be collected at this step.

10 Minute Enzyme Activation Step

The standard protocol for Solaris qPCR Gene Expression Assays incorporates a 15 minute activation step at 95°C. Solaris Assays also perform well using a 10 minute activation step. We recommend using the protocol below to assess the performance of your assay with a shorter activation step.

	Temp.	Time	Number of cycles
Enzyme activation	95°C	10 min	1 cycle
Denaturation	95°C	15 sec	40 cycles
Annealing/Extension	60°C	60 sec	

Fast qPCR Protocol Suggestions

Solaris assays have been developed to give optimal performance under standard thermal cycling conditions. Typically, Solaris assays will also perform well when using a fast qPCR protocol, giving repeatable, specific and sensitive assay results. However, as with any predesigned assay, we recommend that if not following our standard thermal cycling conditions that you check that the performance of your assay (e.g. sensitivity) is not affected. A faster protocol that has worked well in our laboratories, when using Solaris qPCR Gene Expression Master Mix is suggested for your reference below.

	Temp.	Time	Number of cycles
Enzyme activation	95°C	10 min	1 cycle
Denaturation	95°C	5 sec	40 cycles
Annealing/Extension	60°C	20 sec	

Shipping and Storage Conditions

Stability [tests](#) show that Solaris qPCR Gene Expression Assays are stable for several weeks at room temperature. Solaris Assays are shipped on wet ice to provide additional assurance during transit.

For long-term storage we recommend that Solaris qPCR Assays and Master Mixes are stored at -20°C. We also suggest aliquoting into convenient working stock volumes to avoid repeated freeze-thaws of the product.

Following dilution of Solaris qPCR reagents to a final 1X working concentration we recommend running your experiment immediately. However, data generated in our laboratories demonstrate that assays are stable at a 1X concentration for at least 24 hours, providing you with confidence that assay performance will be unaffected within this time-frame (please refer to our technical note entitled "[Stability in the Laboratory](#)").

Ordering Information

AB-4350/INT	Solaris qPCR Gene Expression Master Mix Plus ROX Vial (Intro Pack Size)	100 x 25 µL rxns
AB-4350/A	Solaris qPCR Gene Expression Master Mix Plus ROX Vial	200 x 25 µL rxns
AB-4350/B	Solaris qPCR Gene Expression Master Mix Plus ROX Vial	400 x 25 µL rxns
AB-4350/C	Solaris qPCR Gene Expression Master Mix Plus ROX Vial	1000 x 25 µL rxns

Related Products

Cat No.	Description	Quantity
AX-XXXXXX-XX-XXXX*	Solaris qPCR Gene Expression Assays	100 rxns (1 x 125 µL) 200 rxns (2 x 125 µL) 400 rxns (4 x 125 µL) or 1000 rxns (10 x 125 µL)
AB-1453/A AB-1453/B	Verso cDNA Synthesis Kit	40 x 20 µL rxns 100 x 20 µL rxns
AB-0800/W	ABgene 96-Well PCR Plate (Skirted, Low Profile, White) **	25 plates
AB-2150/W	ABgene Diamond Ultra 384-Well PCR Plate, White **	50 plates
AB-1170	ABsolute™ QPCR Adhesive Seal	50 sheets

* Catalog number is product-specific. Please refer to www.thermo.com/solaris to select the appropriate assay.

** For cyclers compatibility and color choices, see our latest catalog or visit www.thermo.com/pcrplastics

Troubleshooting

For technical information or troubleshooting contact Thermo Scientific Genomics Tech Support:

	<i>Email</i>	<i>Phone</i>
<i>North America (US, Canada, Central/South America)</i>	Techservice.genomics@thermofisher.com	+1 (800) 235-9880
<i>Europe (EU, Middle East, Africa)</i>	Techservice.emea.genomics@thermofisher.com	(+44) 1372 840410
<i>Other Countries</i>	www.thermo.com/dharmacondistributors	

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 6,127,155, 5,677,152 (claims 1 to 23 only) and 5,773,258 (claims 1 and 6 only), and claims outside the US of unexpired patents corresponding to US 4,889,818 and 5,079,352. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim (such as the patented 5' Nuclease Process claims in unexpired counterparts of US Patents Nos. 5,210,015 and 5,487,972 outside of the US) and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA

Use of this product in a passive reference method is covered by the following U.S. Patent: 5,928,907 (claim numbers 12-24, 27-28) and corresponding patent claims outside the US. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product in a passive reference method for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Revised January, 2011. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.
Literature Code: 00024109J01U